

A New York Blood Center Enterprise

EXPANDING OUR ORGANIZATION TO MEET CLINICAL, CELLULAR AND TRANSFUSION PRODUCT AND SERVICE NEEDS FOR PATIENTS. NOW PROVIDING ALMOST ONE MILLION BLOOD PRODUCTS, OVER 450,000 LABORATORY AND MULTI-ASSAY INFECTIOUS DISEASE TESTS AND OVER 12,500 SPECIALTY CLINICAL PROCEDURES ANNUALLY TO HOSPITALE NATIONWIDE.

	COMPREHENSIVE CELL SOLUTIONS	NATIONAL CENTER FOR BLOOD GROUP GENOMICS	NATIO COI BLOO PROGE	NAL DD DD DD
New York Blood Center	Community Blood Center	RENOVATIVE FLOOD RESOURCES	Rhode Island Blood Center	Blood Bank

Objectives

4

- Discuss characteristics of the ABO blood group system, including antigens, antibodies, and genetics.
- Identify and describe several causes of ABO discrepancies, and list techniques used to resolve ABO discrepancies.
- Arrive at appropriate ABO interpretations based on laboratory results.
- Discuss causes of RhD discrepancies, and describe testing used to arrive at appropriate RhD interpretations.



ABO

Case Studies

Objectives

.

- Discuss characteristics of the ABO blood group system, including antigens, antibodies, and genetics.
- Identify and describe several causes of ABO discrepancies, and list techniques used to resolve ABO discrepancies.
- Arrive at appropriate ABO interpretations based on laboratory results.
- Discuss causes of RhD discrepancies, and describe testing used to arrive at appropriate RhD interpretations.



ABO

Case

ABO: Why 1st Blood System Described?

Unique characteristics

- Landsteiner's rule: "naturally occurring antibodies to antigens missing from cells
- Complete agglutinins
- Antibodies reactive over wide thermal range
 RT

Antigen present in large numbers on cells

- Detected without modern technology
 - Centrifugation
 - Antiglobulin reagents

Landsteiner Noble Prize 1933



niversary-of-the-isolation-ofrus/landsteine/2/ Blood Cente

ABO Antigens: the basics

- Levels increase from birth to 3 years of age and then remain stable
- Expressed on RBC, platelets, T and B cells, most body tissues
- Soluble forms in plasma, saliva, most body fluids (depending on secretor status of individual)



ABO Antigen Density

Common blood groups A₁, A₂, B, O and AB • ~ 10⁶/cell = Easily detected by agglutination tests

Adult cells > cord

- -50% fewer antigens on cord cells
 Cord lack branched sugars carrying ABH antigens
 Still easily detected with modern reagents

$A_1 > A_2$

- A₁ five times more antigens than A₂
 Antigens on A₂ cells still easily detected

Weaker subgroups of A or B (AB)

- Rare
 700 50,000/cell
- May not be detectable in agglutination tests
- Cause of ABO discrepancies
 Continuum from strong to weak expression

Community Blood Center

ABO Antibodies

- Appear at 3-6 months; only 15% do not have detectable antibody by 6 months
- · Landsteiner noted "fetal blood contained no agglutinins"
- Titers peak at 5-10 years
- "Naturally occurring" due to environmental stimuli
- Most IgM (except anti-A,B in group O)

Community Blood Center

ABO Genetics

- Over 200 alleles at ABO locus
- Handful convey common blood groups A₁, A₂, B, O, AB in most individuals
 - >99% group A are A₁or A₂
- · Most are rare and associated with weak or null phenotypes
 - · Weak: A subgroups, B subgroups
 - Null: O



ABO Discrepancies

Weak/missing reactivity	•	Extra reactivity
reading reading		Extra reading

Front type
Back type

Front type
 Back type

Daonty			Duontypo	
	From	t type	Back	к Туре
	Anti-A	Anti-B	A1 cell	B cell
Weak	1+	0	0	4+
reactivity	0	0	2+	4+
Missing	0	0	0	4+
reactivity	0	4+	0	0
Additional	1+	1+	4+	4+
reactivity	4+	0	2+	4+

Community Blood Center

ABO Discrepancies: Where to begin?

Correct sample?? Technical error??

Patient History

- · Very young/old: weak ABO antibodies
- · Alloantibody that might interfere with reverse grouping
- Strong cold autoantibody
- May interfere with both forward and reverse grouping
- Bone Marrow Transplant
- Recent transfusion
- Diagnosis
 - · Weak antigens in leukemia, pregnancy, cord samples



Objectives

- Discuss characteristics of the ABO blood group system, including antigens, antibodies, and genetics.
- Identify and describe several causes of ABO discrepancies, and list techniques used to resolve ABO discrepancies.
 Arrive at appropriate ABO interpretations based on laboratory results.
- Discuss causes of RhD discrepancies, and describe testing used to arrive at appropriate RhD interpretations.

RhD discussion & Cases
& Cases
Community

ABO

Case

Studies

Case #1









What's anti-A1_lectin?





What's the difference between A_1 and A_2 phenotypes?

• Frequency in European population

- ~80% group A individuals A₁
- ~20% group A individuals A₂
- Other A subgroups rare (A₃, A_x, A_{el}, etc.)

Antigen differences

- · Quantity of A antigens on cells
 - A₁ cells have approximately 5 times as many A antigens as A₂ cells
- Qualitative differences
 - Antigens of A1 individuals more branched
 - Why A subgroup individuals can make anti-A1.

Fung MK, Eder AF, Spitalnic SL, Westhoff CM. Technical Manual. 19th ed. Bethesda, MD: AABB; 2017: 271-272

Community Blood Center

More about anti-A1

- Detected in 1-8% of A₂ individuals, 22-35% of A₂B individuals
 Often "naturally occurring"
- Usually IgM, reacts best at room temperature or below.
 - Generally not considered clinically significant
- Reports in literature of hemolytic anti-A1
- Helmich F, et al. Acute hemolytic transfusion reaction due to a warm reactive arti-A1. Transfus
 Transfusion recommendations: XM compatible units

Fung MK, Eder AF, Spitalnic SL, Westhoff CM. Technical Manual. 19th ed. Bethesda, MD: AABB; 2017: 274

- A₂ RBCs
- O RBCs

Community Blood Center

2018;58;1163-1170.





Does the patient's plasma contain cold autoantibody?

				Rh			Ke	:II	Du	ıffy	Ki	dd		M	NS			Results	
		D	C	Ε	с	e	К	k	Fv³	Fv ^b	Jka	Jkb	M	N	S	s	5' RT	LISS	LISS
			1	nc	nr	02	acti	ive	<u>م</u>	الص	at	RT	5 I					37C	IAT
1		+		110		00	101		- 0	Cii	а		•	+	+	+	1+	0	0
2		+			1-	ne	eqa	ati	ve	ce	11			+	0	+	1+	0	0
3		+	0 + + 1+												1+	0	0		
4		+	0	0	+	+	0	+	2	2	+	0	+	+	0	+	1+	0	0
5		0	+	0	+	+	0	+	+	0	×	0	+	+	0	0	1+	0	0
6		0	0	+	+	+	0	+	0	+	+	7	0	+	0	+	1+	0	0
7		0	0	0	+	+	+	+	0	+	+	0	1	0	+	+	1+	0	0
8		0	0	0	+	+	0	+	+	+	0	+	0	+	+	+	1+	0	0
9		0	0	0	+	+	0	+	+	+	0	+	+	0	0	Ż	1+	0	0
10	1-	+	+	0	0	+	0	+	+	0	+	+	+	+	+	0	0	0	0
11		+	0	0	+	+	+	+	0	0	+	+	0	+	+	+	1+	0	0
Auto																	1+	0	0
	Autocontrol reactive at 5' RT													munity					

COLD SCREEN

	Cold Screen	
	30' RT	30' 4C
SCI	2+	4+
SCII	2+	4+
I-negative	0	3+
Auto	2+	4+

What is pre-warming?

- Incubate all reagent red cells and patient plasma (and pipet) at 37C
 prior to testing
- Add patient plasma to reagent cells quickly without changing the temperature of the testing environment from 37C
- Incubate all tubes at 37C (~30 min)
- DO NOT CENTRIFUGE! (centrifugation will quickly cool the sample in the tubes)
- · Shake and read settled tubes



One word of caution!

NEVER utilize pre-warm testing unless you know what you are pre-warming!!

(demonstrate that the patient has cold autoantibody prior to prewarming)





Scenario #3

	Front	Туре			E	Back Typ	е	
Anti-A	Anti-B	Anti-D	Rh Control	A ₁	A ₂	В	0	Auto control
4+	0	4+	0	1+	1+	4+	1+	1+

These same results (same as scenario #2) may be due to rouleaux!!

- In vitro phenomenon due to abnormal patient plasma protein concentration
- Can be seen in any test involving patient plasma (including back type)
 "Stack of coins," refractile aggregation of RBCs
 Can be mistaken for agglutination macroscopically

Fung MK, Eder AF, Soltahic SL, Westhoff CM, Technical Manual, 19th ed. Bethesda, MD: AABB: 2017: 370-371

Community Blood Center

Can you see the difference?



Steps to saline replacement...





Steps to saline replacement...











Cells for back type: Pooled





Identifying the alloantibody

			Rh			К	ell	Du	ıffy	Ki	dd		М	NS			Result	s
	D	С	E	с	e	к	k	Fya	Fyb	Jka	Jkb	м	N	S	s	5' RT	LISS 37C	LISS IAT
1	+	+	0	0	+	0	+	+	+	+	+	+	+	+	+	1+	0	0
2	+	+	0	0	+	+	+	0	+	0	+	0	+	0	+	0	0	0
3	+	0	+	+	0	0	+	+	0	+	+	+	0	+	+	2+	1+	1+
4	+	0	0	+	+	0	+	0	0	+	0	+	+	0	+	1+	0	0
5	0	+	0	+	+	0	+	+	0	+	0	+	+	0	0	1+	0	0
6	0	0	+	+	+	0	+	0	+	+	+	0	+	0	+	0	0	0
7	0	0	0	+	+	+	+	0	+	+	0	+	0	+	+	2+	1+	1+
8	0	0	0	+	+	0	+	+	+	0	+	0	+	+	+	0	0	0
9	0	0	0	+	+	0	+	+	+	0	+	+	0	0	+	2+	1+	1+
10	+	+	0	0	+	0	+	+	0	+	+	+	+	+	0	1+	0	0
11	+	0	0	+	+	+	+	0	0	+	+	0	+	+	+	0	0	0
Aut o																0	0	0

Scenario #4

	Front	Туре		Back Type						
Anti-A	Anti-B	Anti-D	Rh Control	A ₁	A ₂	В	0	Auto control		
4+	0	4+	0	1+	1+	4+	1+	0		

Once the antibody is identified, resolve the typing discrepancy by...

- · Prewarming the back type
- · Using RBCs for your back type testing that don't express the
 - For example, M-negative A₁ cells & M-negative B cells
 Enzyme-treated cells

Community Blood Center

Review of the 4 scenarios:





Review of the 4 scenarios:





Review of the 4 scenarios:



Ca	ase	#2						
	Front	Туре			E	Back Typ	e	
Anti-A	Anti-B	Anti-D	Rh Control	A ₁	A ₂	В	0	Auto control
0	0	4+	0	0	0	0	0	0

Ways to promote/strengthen reactivity of back type:

- Increase incubation time at 22C
- · Decrease temperature*
- Increase plasma:cell ratio (use 4 drops of plasma & 1 drop of cells in each tube)

* Be careful: Many individuals have cold autoantibodies!

Community Blood Center



Wait... What is the Rh control?



Strong cold agglutinins



Community Blood Center





Mixed field



Case #4

Patient is group A+ and received O- blood Patient is group A+ and received A- blood

	Front	Туре			E	Back Type	e				
Anti-A	Anti-B	Anti-D	Rh Control	A ₁	A ₂	В	0	Auto control			
4+ ^{mf}	0	4+	0	0	0	3+	0	0			
What	What is the best explanation for these results?										
Patier	Patient is group O+ and received A+ blood										
Patie	Patient is group A+ and received O+ blood										

Community Blood Center

Ca	ase	#4	k i					
	Front	Туре			E	Back Typ	е	
Anti-A	Anti-B	Anti-D	Rh Control	A ₁	A ₂	В	0	Auto control
4+ ^{mf}	0	4+	0	0	0	3+	0	0
				L	L	ooks l aroup	ike A	









Case #4

	Front	Туре			E	Back Typ	е	
Anti-A	Anti-B	Anti-D	Rh Control	A ₁	A ₂	В	0	Auto control
4+ ^{mf}	0	4+	0	0	0	3+	0	0

Patient is group A+ and received O+ blood

Important things to know!!!

- Transfusion of non ABO-identical RBCs affects the <u>FRONT</u> type
 Transfusion does not usually interfere with the back type
 Use the mixed-field reactions to determine what type of blood patient
- received

Community Blood Center

Weak A subgroup



What do we mean by adsorption/elution?





Weak A subgroup RBCs

Community Blood Center

What do we mean by adsorption/elution?



Recent massive transfusion with group **O RBCs**

In cases of massive transfusion, it is possible to entirely replace patient blood volume with transfused cells



These results could indicate a group A patient (according to back type) who
recently has been massively transfused with group O RBCs

Get patient history!!!

Hematopoietic Stem Cell/ **Bone Marrow Transplant**

Front Type					E	Back Typ	е	
Anti-A	Anti-B	Anti-D	Rh Control	A ₁	A ₂	В	0	Auto control
0	0	3+	0	0	0	4+	0	0

- Patient's blood type may change
 Example: Group A patient receives transplant from group O donor
 May receive patient sample during engraftment
 Both donor and recipient cell populations detected
 Mixed field results
 It transplant encourses the security

 - If transplant successfully engrafts
 Patient will assume donor blood type
 May not develop antibodies to antigens his/her cells formerly expressed

Community Blood Center

	Revi	iew: Common ABC) Discrepancies
	Discrepancy	Cause	Resolution
	Extra reactions	A subgroup with anti-A1	Test plasma with A ₂ cells, test cells with anti-A1 lectin
		Rouleaux	Confirm rouleaux microscopically, saline replacement
ď		Cold autoantibody	Pre-warm back type
ack ty		Cold alloantibody	Identify the alloantibody, then pre-warm back type or use antigen negative cells for back type
ä	Missing reactions	Neonate	Back type not performed <4 months
		Immunosuppression	Incubate back type at room temperature, increase plasma:cell ratio, decrease temperature
	Extra reactions	Cold agglutinin coating red cells	Wash cells with warm saline
type	Missing reactions	Weak A or B subgroups	Test cells with anti-A,B, adsorption/elutions, genomic testing
Ŧ	Mixed-field	Multiple cell populations in	Acquire patient history
-o	reactions	sample due to recent	
÷		transfusion/Bone marrow	
		transplant	
	Special	Bone Marrow Transplant	Acquire patient history
cor	nsiderations	 Recent massive transfusion with group O RBCs 	

Review: ABO discrepancies

	[Dossik	lo Evo	lanati	ane			M	lavs to	Rosol	ve	
		Possib Cold Au Recent ¹ Cold rea A subgro BMT (gr BOUSOL Immuno	I e Exp toantibo Transfus acting all oup with oup A+	lanatio dy ion oantibo n anti-A1 to O+) ision	dy			In BC Ni Sa W Di Id Te	Vays to crease in guire ab o further line rep arm was acrease i entify al est with	Resol ncubatio back typ out pati testing lacemen in front incubatio loantibo anti-A1 l	ve n time of back type e ent history required t of backtype, type n temperature dy ectin	•
1	Front	Type:	Test	ng pa	tient	В	ack T	ype:	Testii	ng		
	Front	Type:	Testi cells	ng pa	tient	В	ack T	ype: ent pla	Testii	ng		
	Front	Type:	Testi cells Anti- A,B	ng pa	Rh Cont.	A1 cells	ack T patie	ype: ent pla	Testii asma o cells	Auto Cont.	Possible Explanation	Ways to Resolve
1	Front Anti- A 4+	Type: Anti- B 0	Testi cells Anti- A,B 4+	ng pa	Rh Cont. 0	A1 cells 1+	ack T patie A2 cells 0	ype: ent pla cells 4+	Testin asma o cells 0	Auto Cont. 0	Possible Explanation	Ways to Resolve
1	Anti- A 4+	Type: Anti-B 0 4+	Testi cells Anti- A,B 4+ 4+	ng pa	Rh Cont. 0	A1 cells 1+ 4+	ack T patie ^{A2} cells 0 3+	ype: ent pl: cells 4+	Testin asma o cells 0	Auto Cont. 0	Possible Explanation •	Ways to Resolve •



Objectives

.

•

- Discuss characteristics of the ABO blood group system, including antigens, antibodies, and genetics.
 - Identify and describe several causes of ABO discrepancies, and list techniques used to resolve ABO discrepancies. ABO Case Studies based on laboratory results.
- Discuss causes of RhD discrepancies, and describe testing used to arrive at appropriate RhD interpretations.



RhD discrepancies

Causes:

- Incorrect sample/technical error
- Unusual RhD antigen expression
 - Weak D
- Partial D
- Difference in antibody clones in different reagents
- Antibodies to different epitopes of D antigen
 Different methodologies
- Exampe: IS vs IAT
- Patient history
- BMT? Transfusion?

Community Blood Center

More on RhD discrepancies

Narrow specificity of monoclonal Rh reagents

Single epitope vs entire antigen

Expected that cells with unusual antigens (RhD variants) will react differently with different monoclonal antisera

- · Strong, weak or negative reactions
- Patients/donors may "change" Rh type if different reagent used

RhD discrepancies: strategies

Patient race

Different RhD variants are more prevalent in ٠ individuals in specific ethnic groups

Test with multiple reagents containing different anti-D clones and human polyclonal anti-D

- · Known reaction pattern with some common partial Ds
- Run controls
- Record lot number and reagent ID (e.g., 114757 IM Series 4)

 - Human polyclonal anti-D
 IM Series 4 and 5 anti-D
 Gammaclone anti-D
 Ortho Bioclone anti-D

Community Blood Cente

RhD discrepancies: strategies

Test with Partial RhD Typing Kit

- · Panel of 12 monoclonal anti-Ds
- Classify 16 different partial or weak D depending on reactivity



Test for Cc, Ee, and low prevalence antigens, if applicable

See FactsBook

Molecular testing

Definitive way to accurately determine RhD status

Community Blood Center

Case 1

Memorial Hospital reports a patient previously typed as Rh positive in 2002 now appears Rh negative

Questions

- Race of patient?
- History of BMT???
- What reagent/method used in 2002?
- · What reagent/method used currently?

Case 1

Memorial Hospital reports a patient previously typed as Rh positive in 2002 now appears Rh negative

Questions

- Race of patient? African American
- History of BMT??? No
- What reagent/method used in 2002? Gammaclone anti-D/tube-immediate spin reading only
- What reagent/method used currently? Ortho gel

Community Blood Center

IRL Testing

Different Anti-D reagents									
	IM Gamma- clone IS/IAT	IM Series 4	IM Series 5	Ortho Bioclone	Human Poly- clonal IS/IAT				
Patient	Pos/Neg	Neg	Neg	Neg	Neg/Neg				
Pos cont	Pos	Pos	Pos	Pos	Pos				
Neg cont	Neg	Neg	Neg	Neg	Neg				

Ortho Gel not available

(Other Rh Antigens					
	Anti-C	Anti-E	Anti-c	Anti-e		
Patient cells	0	0	4+	4+		

Reactivity of FDA licensed Anti-D with Unusual (Variant) D Antigens

Reagent	lgM	lgG	Normal D Pos IS test	DVI IS/weak D	DBT IS/weak D	D ^{Har} Whites IS/weak D	Crawford Blacks IS/weak D
Gammaclone	GAMA401	F8D8	Pos	Neg/Pos	Pos	Pos	Pos /Neg
Immucor Series 4	MS201	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Immucor Series 5	Th 28	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Ortho Bioclone	MAD2	Poly- clonal	Pos	Neg/Pos	Neg/Pos	Neg/Neg	Neg
Ortho Gel	MS201	NA	Pos	Neg	Pos	Pos	Neg
Polyclonal			Pos	Neg/Pos	Neg/Neg	Neg/Neg	Neg/Neg



Reactivity of FDA licensed A same pattern with Unusual (Variant) D Antiger

Reagent	IgM	lgG	Normal D Pos IS test	DVI IS/weak D	DBT IS/weak D	D ^{Har} Whites IS/weak D	Crawford Blacks IS/weak D
Gammaclone	GAMA401	F8D8	Pos	Neg/Pos	Pos	Pos	Pos /Neg
Immucor Series 4	MS201	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Immucor Series 5	Th 28	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Ortho Bioclone	MAD2	Poly- clonal	Pos	Neg/Pos	Neg/Pos	Neg/Neg	Neg
Ortho Gel	MS201	NA	Pos	Neg	Pos	Pos	Neg
Polyclonal			Pos	Neg/Pos	Neg/Neg	Neg/Neg	Neg/Neg

Community Blood Center

ANTIGEN

FactsBook

http

Case 1: Review

African American patient

Reactive with Gammaclone anti-D

- IgM clone strongly reactive at IS with low incidence Crawford, not D antigen, on cells
 IgG clone negative, i.e., weak D test neg
- Nonreactive with Ortho Gel

Other reagents tested consistent with Crawford Crawford seen with Ce/ce and ce/ce phenotypes Consider Rh negative for transfusion

 Can make multiple Rh system antibodies – see FactsBook

We would confirm this with *RHD* genotypeuse serology and molecular testing together!

Case 2

A donor's historical ABO/Rh is A negative from 3 donations in the 1990s. She was tested at her doctor's office and told she was Rh positive

Questions

- Race of Donor?
- History of BMT???
- What reagent/method used in 1990s?
- · What reagent/method used currently?

Case 2

A donor's historical ABO/Rh is A negative from 3 donations in the 1990s. She was tested at her doctor's office and told she was Rh positive

Questions

- Race of Donor? Caucasian
- History of BMT??? No
- What reagent/method used in 1990s? Immucor human polyclonal anti-D/PK7300
- What reagent/method used currently? Ortho gel

Community Blood Center

IRL Testing

Different Anti-D Reagents									
	IM Gamma- clone	IM Series 4	IM Series 5	Ortho Bio- clone IS/IAT	Human Poly- clonal IS/IAT				
Patient	Pos	Pos	Pos	Neg/Neg	Neg/Neg				
Pos cont	Pos	Pos	Pos	Pos	Pos				
Neg cont	Neg	Neg	Neg	Neg	Neg				

Ortho Gel not available

Anti-C Anti-E Anti-c Anti-e Patient 0 0 4+ 2+ cells	Other Rh Antigens									
Patient 0 0 4+ 2+ cells		Anti-C	Anti-E	Anti-c	Anti-e					
	Patient cells	0	0	4+	2+					

Community Blood Center

Reactivity of FDA licensed Anti-D with Unusual (Variant) D Antigens

Reagent	lgM	lgG	Normal D Pos IS test	DVI IS/weak D	DBT IS/weak D	D ^{Har} Whites IS/weak D	Crawford Blacks IS/weak D
Gammaclone	GAMA401	F8D8	Pos	Neg/Pos	Neg/Pos Pos		Pos /Neg
Immucor Series 4	MS201	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Immucor Series 5	Th 28	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Ortho Bioclone	MAD2	Poly- clonal	Pos	Neg/Pos	Neg/Pos	Neg/Neg	Neg
Ortho Gel	MS201	NA	Pos	Neg	Pos	Pos	Neg
Polyclonal			Pos	Neg/Pos	Neg/Neg	Neg/Neg	Neg/Neg

Reactivity of FDA license Same pattern with Unusual (Variant) D Art A

Reagent	lgM	lgG	Normal D Pos IS test	DVI IS/weak D	DBT IS/weak D	D ^{Har} Whites IS/weak D	Crawford Blacks IS/weak D
Gammaclone	GAMA401	F8D8	Pos	Neg/Pos	eg/Pos Pos Pos		Pos /Neg
Immucor Series 4	MS201	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Immucor Series 5	Th 28	MS26	Pos	Neg/Pos	Pos	Pos	Neg
Ortho Bioclone	MAD2	Poly- clonal	Pos	Neg/Pos	Neg/Pos	Neg/Neg	Neg
Ortho Gel	MS201	NA	Pos	Neg	Pos	Pos	Neg
Polyclonal			Pos	Neg/Pos	Neg/Neg	Neg/Neg	Neg/Neg

Case 2: Partial D Typing Kit Results

KelD	Anti-D Cell Line	Weak D Type 1 and 2 ⁸	DI & DNU	DII	DIV	٥v	DCS	DM	DVII	DOL	OFR	DNH	DAR	DAR-E	DHK ⁴ 8 Dau-4	DBT	Ro ^{Rer} l	Pos Cont.	Neg Cont	Case 2
A	LHM76/58	+	+	٠	+	+10	+	0	+	+	٠	٠	+	0	0	0	(+)0	3+	0	0
В	LHM76/59	٠	+	٠	0	٠	+	٠	+	+	٠	٠	٠	٠	+	0	0	3+	0	0
С	LHM174/102	(+)0	+	٠	0		+		+	0	0	٠		0	0	0	0	3+	0	0
D	LHM50/28	٠	+	٠	+	٠	+	0	+	+	٠	٠	٠	٠	+	0	0	3+	0	0
Ε	LHM169/81	•	+	٠	0	0	+	0	+	+	٠	٠		0	0	0	0	3+	0	0
F	ESD1	٠	+	٠	0	٠	+	•	+	+	٠	٠	٠	٠	+	0	0	3+	0	0
G	LHM76/55	•	+	٠	0	٠	+	+	+	+	٠	٠	•	٠	•	0	0	3+	0	0
н	LHM77/64	٠	0	٠	0	٠	+	•	+	+	٠	٠	٠	٠	-0	0	0	3+	0	0
1	LHM7045	(+)0	+	٠	0	0	0	0	+	0	0	0		0	0	0	0	3+	0	0
J	LHM59/19	•	+	٠	+	+	+	0	0	0	0	(+)		(+)	+	٠	0	3+	0	0
K	LHM169/80	٠	+	٠	+	٠	+	0	+	+	٠	٠	٠	٠	0	0	0	3+	0	0
L	LHM57/17	+	+	+	+	+	0	0	+	+	0	•	•	0	0	+	0	3+	0	0

Community Blood Center

Community Blood Center

Case 2: Review

Caucasian patient/donor

Testing consistent with RH33/DHar

- Reactive with IM Gammaclone anti-D, Series 4, Series 5, Ortho Gel
- Nonreactive with Human polyclonal
- Partial D Typing Kit RhD^{Har}

RH33/D^{Har} seen with c+, weak e+; see *FactsBook* Consider Rh negative for transfusion

Blood center follow-up:

Former donor units labeled as RhD-negative

We would confirm this with *RHD* genotypeuse serology and molecular testing together!



amazon.com/Bio 3p/0124158498

RhD discrepancies: take home message

- · Most donor/patients carry common Rh genes, express normal antigen on cells and react well with all commercial anti-D reagents
- Cells from individuals with unusual genes, expressing unusual antigens will react differently with different monoclonal antisera
- Serologic testing and molecular testing used together to accurately determine RhD status •



Community Blood Center

Objectives

- Discuss characteristics of the ABO blood Introduction group system, including antigens, antibodies, and genetics.
- . Identify and describe several causes of ABO discrepancies, and list techniques used to resolve ABO discrepancies.
- Arrive at appropriate ABO interpretations based on laboratory results.
- Discuss causes of RhD discrepancies, and describe testing used to arrive at appropriate RhD interpretations.



ABO

Case Studies

Community Blood Cente

References

Fung MK, Eder AF, Spitalnic SL, Westhoff CM. Technical Manual. 19th ed. Bethesda, MD: AABB; 2017: 274.

Helmich F, et al. Acute hemolytic transfusion reaction due to a warm reactive anti-A1. Transfusion. 2018;58;1163-1170.



A New York Blood Center



Blood Bank

SAVEALIFENOW.ORG